STRETTA OVERVIEW

The Stretta System consists of the Stretta Catheter — a single-use flexible device with needle electrodes, and the Mederi RF Generator. The Stretta Procedure is typically performed in an endoscopy unit or ambulatory surgery. In the hands of a physician Stretta is used to treat the muscle of the lower esophageal sphincter (LES) and gastric cardia in patients with Gastroesophageal Reflux Disease (GERD). The Stretta Catheter is positioned at the junction of the esophagus and stomach (GE junction). The integrated balloon is inflated and needle electrodes are then extended from the device to deliver radiofrequency energy into the target tissue.

Federal law restricts this device to sale by or on the order of a physician.

MEDERI THERAPEUTICS INC

stretta-therapy.com

800 Connecticut Ave.
Norwalk, CT 06854
203-930-9900

MediMark Europe
11, rue Emile Zola, BP 2332
38033 Grenoble Cedex 2, France

174-5575G

English

INDICATIONS FOR USE

The Stretta System is intended for general use in electrosurgical coagulation of tissue and intended for use specifically in the treatment of gastroesophageal reflux disease (GERD).

CONTRAINDICATIONS

There are no known absolute contraindications to the use of radiofrequency in humans. The use of the Stretta System is contraindicated when, in the judgment of the physician, radiofrequency surgical procedures would be contrary to the best interests of the patient.

The following is a list of patient groups in whom the use of the Stretta System for the treatment of GERD may be contraindicated:

1. Subjects under the age of 18
2. Pregnant women
3. Patients without a diagnosis of GERD
4. Hiatal hernia > 2 cm
5. Achalasia or incomplete LES relaxation in response to swallow
6. Poor surgical candidate, ASA/IV classification

WARNINGS

For the Stretta System, only the Mederi RF Generator can be used with the Stretta Disposable Catheter.

Carefully read all instructions for use for Stretta Catheter and Mederi RF Generator.

As with any procedure, use of sedation or general anesthesia for patients with airway or cardiopulmonary compromise (such as obesity, airway abnormality and sleep apnea) may result in increased risk of complications.

This product must be disposed of according to local laws and regulations.

WARNINGS RELATED TO THE USE OF ELECTROSURGERY

- Use of electrosurgery for patients with internal or external pacemakers or ICDs, neural stimulators, implantable drug delivery systems, or any other electrical implantable devices and monitoring equipment may require special consideration. The appropriate physician and/or manufacturer should be consulted prior to surgery.
- This electrosurgical equipment is intended for use only by qualified medical personnel trained in the use of the Stretta System.
- Burns to the physician's hands are possible if the electrode comes into contact with a metal instrument or surface.
- Interference produced by the operation of high-frequency surgical equipment may adversely influence the operation of other electronic medical equipment such as monitors and imaging systems.
- Do not stack equipment on top of the Mederi MDRF1 Generator or place the generator on top of or adjacent to electrical equipment. This may block access to the unit and not allow for proper ventilation.
- Provide as much distance as possible between the Mederi MDRF1 Generator and other electronic equipment (such as monitors). An activated generator may cause interference with them.
- If it is necessary to locate the Mederi MDRF1 Generator adjacent to or stacked with other equipment, then observe the generator and other equipment for proper operation.
- Do not position the Generator so that it is difficult to disconnect the power cords.
- Dispose of Stretta disposable Catheter and accessories in an environmentally safe manner according to local regulations.
- CAUTION: Electrical shock hazard. Do not remove the cover of the Mederi RF Generator. Refer servicing to Mederi Therapeutics qualified personnel. There are no user serviceable parts inside the Generator.
- CAUTION: No unauthorized modification to this device is allowed.
- CAUTION: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- Do not perform the procedure if flammable or explosive media are present, such as flammable anesthetics, skin preparation agents, or bio-intestinal gases.
- DANGER: Risk of explosion if used in the presence of flammable anesthetics.
- WARNING: The risk of igniting flammable gases or other materials is inherent in the application of RF power. These substances may come in the form of anesthetics, oxygen, cleaning or disinfection agents, adhesive solvents, endogenous gases, surgical drapes, etc. and should be removed from the area of RF energy delivery. Non-flammable agents should be used for cleaning and disinfection wherever possible as a means to reduce the possibility of ignition. Flammable agents used for cleaning or disinfecting, or as solvents of adhesives, should be allowed to evaporate before the application of RF power to reduce the chance of ignition.
- WARNING: Flammable agents used for cleaning or disinfecting, or as solvents of adhesives, should be allowed to evaporate before the application of RF surgery.
- WARNING: There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before RF surgical equipment is used. Attention should be called to the danger of ignition of endogenous gases. Some materials for example cotton, wool and gauze, when saturated with oxygen may be ignited by sparks produced in normal use of RF surgical equipment.
- Stimulation of nerves and muscles. A known risk of high-frequency surgery is the unintentional electric stimulation of the patient's nerves and muscles. This stimulation can result from low-frequency electrical currents that are caused either by low-frequency current sources or due to electrical arcs between an active electrode and the patient's tissue. Electric alternating current with a frequency above 300 kHz is unable to stimulate nerves and muscles. Note that the electric alternating current used by this generator is 460kHz. During controlled application of radiofrequency to muscles in the GI tract, the unavoidable electric arcs between an active electrode and the tissue nevertheless have the effect that

a portion of the high-frequency alternating current is rectified, from which more or less strongly modulated, low-frequency current components result which stimulate electrically stimulative structures such as nerves and muscles.

This can result in more or less strong spasms or muscle contractions.

WARNING: When using high-frequency energy on electrically stimulative structures, contractions of the affected muscles must be taken into account. This can occur, for example, during endoscopic operations in the area of the lower esophageal sphincter or the internal anal sphincter muscles.

WARNINGS RELATED TO THIS SPECIFIC PROCEDURE

- Inspect all components prior to use for any obvious signs of damage that may have occurred during transit and/or storage. In particular, the electrode cable should be checked prior to use.
- Ensure adequate volume of irrigation fluid is available prior to start of procedure.
- Do not attach suction tube until the electrode catheter is positioned in the tissue.
- Be sure that the reusable cable has been cleaned prior to use.
- Do not attempt to re-sterilize and/or reuse the electrode catheter.
- Reuse of Mederi single-use items is contraindicated in all patients. Like all single-use items, use on more than one patient may result in diseases of cross-contamination and possible transmission of blood borne pathogens.
- Balloon must be in the deflated state during insertion, treatment location adjustment, and removal of the catheter from target organs.
- Pressure relief valve should be used for antegrade lesions only. Do not use for pull-back treatments.
- If pressure relief valve malfunctions in use, remove it and replace.
- Needle electrodes must be in the retracted position during insertion, treatment location adjustment, and removal of the catheter from target organs.
- WARNING: To avoid overlap of treatment sites ensure that treatment sites are spaced a minimum of 5mm apart in all directions.
- WARNING: Failure to follow proper shutdown sequence could cause damage to the Generator.

CAUTION: Use of accessories not outlined in this Operator's Manual could affect the safety of the Stretta System. Use only power cord supplied by Mederi.

CAUTION: "Spring-tipped" guidewires are not recommended.

CAUTION: Due to elastic properties of the Irrigation Tubing material, excessive longitudinal force may result in stretching (flat tube on the pump rotor) and as a result, reduction of the irrigation flow.

CAUTION: Failure to achieve good skin contact by the entire adhesive surface of the return electrode may result in a burn or poor electrical performance. Excessive hair, moist skin, lotions, creams or skin creases, may result in inadequate electrical contact.

CAUTION: Do not use steam or heat sterilization on the Mederi Generator. Do not soak in glutaraldehyde or other disinfectant. Do not allow any liquid to pass into any electrical connections or into the interior of the Generator.

When conducting periodic inspections of the Mederi Generator, refer to the chapter, Care & Maintenance, if testing reveals a defect that could harm the patient, employees, or third parties, immediately notify your Mederi representative. The device should not be used until it has been properly repaired/serviced by an authorized Mederi technician or representative.

The following is a list of patient groups who have not been specifically studied using the Stretta System for the treatment of GERD. The performance characteristics of the device for the following groups of patients have not been established.

- Patient has an implant near the LES that could be conductive with RF energy
- Normal 24-hour pH study
- ALL GERD symptoms completely unresponsive to properly dose-escalated anti-secretory medication
- Barrett's metaplasia
- Poor surgical candidate
- Presence of dysphagia, esophageal bleeding, or gas bloat
- Active esophagitis grades III or IV by Savary criteria
- Endocarditis risk (mitral valve prolapse, heart valve replacement, etc.)
- Untreated or unstable hypertension, diabetes mellitus, heart disease, collagen vascular disease, steroid use, immunosuppressed state, or cardiac pacemaker
- Abnormal blood coagulation or use of anticoagulant or platelet anti-aggregation therapy.

These complications are rarely seen but could potentially occur with the use of electrosurgery for the treatment of GERD (alphabetical order).

- Bleeding - transient
 - Bloating
 - Chest pain - transient
 - Difficulty belching - transient
 - Dysphagia - transient
 - Epigastric discomfort - transient
 - Esophageal mucosal laceration
 - Fever - transient
 - Injury to esophageal mucosa
 - Perforation
 - Pharyngitis
 - Vomiting - transient with potential for bleeding or Esophageal injury. If any vomiting occurs, contact your treating physician immediately. Excessive vomiting may result in perforation and more serious injury resulting in death.
- The following complications have not been seen, but could possibly occur infrequently (alphabetical order).
- Achalasia

- Delayed gastric emptying - transient
- Dental injury
- Dyspnea
- Infection
- Larynx injury
- Worsened GERD

WARNING: Use only sterile, chilled, plain water as irrigation fluid for this device. Do not use saline or other conductive/ionic/volatile solution, as this may lead to improper device performance. Do not expand balloon to more than 15 psi or 30 ml of air.

PRECAUTIONS RELATED TO THE USE OF ELECTROSURGERY

- Safe use of monopolar electrosurgery demands proper grounding of the patient. Follow the manufacturer's directions for the placement of the dispersive electrode and for proper insulation between the patient and any metal surfaces.
- Maintain safe handling techniques when Catheter is in use due to electrical field and potentially hot needle electrodes.
- Do not touch the tip of the needle electrode and the dispersive electrode at the same time - especially while operating the Mederi RF Generator. Superficial skin burns could occur.
- Ensure that the dispersive electrode pad achieves complete contact with cleaned and hair-less patient skin.
- Activation of an electrosurgical device when not in contact with target tissue or in a position to deliver energy to target tissue, but when in close approximation to a metal object, may result in capacitive coupling which could result in shock or discharge.
- Simultaneous use of electrosurgical electrodes with conductive irrigation fluid (such as saline) may alter the path of electrical energy away from target tissue. Do not operate the Stretta Disposable Catheter in the presence of saline or other conductive irrigation fluid.
- This Mederi RF Generator is in compliance with current required standards for electromagnetic interference and should not present problems to other equipment or be affected by other equipment. As a precaution, avoid use of the Mederi RF Generator in close proximity to other equipment.
- Interference from other electrical equipment may cause interference in the operation of the Mederi MDRF1. Also, interference to other electrical equipment may be caused by the Mederi MDRF1. If interference occurs, increase the physical separation distance between the Mederi RF Generator and other electrical equipment. Avoid the overlap of cables from the Mederi MDRF1 Generator with cables from other electrical equipment.
- Power Cords: Do not wrap power cords around metal objects. This may induce currents that could lead to shock, fire, or other injury to the patient or user. All power cords should be positioned in a way to avoid contact with the patient or other cables.
- Emissions produced by the operation of the Mederi RF Generator together with the Stretta Disposable Catheter may adversely influence the operation of other electronic medical equipment such as monitors and imaging systems.
- The immunity of the Mederi RF Generator to interference by other equipment may be decreased when the Stretta Disposable Catheter is used together with the Mederi RF Generator.
- When conducting periodic inspections of the Mederi RF Generator, refer to the Care and Maintenance section in this Operator's Manual. If testing reveals a defect that could harm the patient, employees, or third parties, immediately notify your Mederi representative. The device should not be used until it has been properly repaired/serviced by an authorized Mederi Therapeutics technician.
- Do not operate the RF Generator when the ambient temperature is below 15° C.
- The entire area of the neutral electrode should be reliably attached to the patient's body and as close to the operating field as possible.
- The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth (for example operating table supports, etc.). The use of antistatic sheeting is recommended for this purpose.
- Skin-to-skin contact (for example between the arms and body of the patient) should be avoided, for example by insertion of dry gauze.
- When HF surgical equipment and physiological monitoring equipment are used simultaneously on the same patient, any monitoring electrodes should be placed as far as possible from the surgical electrodes. Needle monitoring electrodes are not recommended. In all cases, monitoring systems incorporating high frequency current limiting devices are recommended.
- The cables to the surgical electrodes should be positioned in such a way that contact with the patient or other leads is avoided. Temporarily unused active electrodes should be stored in a location that is isolated from the patient.
- Non-flammable agents should be used for cleaning and disinfection wherever possible (see Care and Maintenance). Flammable agents used for cleaning or disinfecting, or as solvents of adhesives, should be allowed to evaporate before the application of HF energy.
- Failure of HF surgical equipment could result in an unintended increase of output power.

REQUIRED ITEMS FOR PERFORMING STRETITA THERAPY

1. Mederi RF Generator, accessories and Operator's Manual
2. Compatible return electrode pad (not included, see Operator's Manual)
3. Stretta Catheter and accessories (guidewire, syringe, PRV, irrigation tubing)
4. Chilled bag of sterile water (Water only, DO NOT USE SALINE)
5. Standard Endoscope

PATIENT PREPARATION

The physician may elect to modify these suggested steps for patient preparation:

1. Prepare the patient for standard upper endoscopy, placing the patient in the left lateral decubitus position.
2. Patient should have neck extended with bite block. Tilt bed slightly (head higher) to avoid risk of aspiration.
2. Prepare the patient using standard technique for monopolar electrosurgery.

3. Apply the return electrode pad to a cleaned and hairless area on the patient's right mid scapular area of the mid line. Area must be clean and hairless to ensure proper electrical contact.
4. Patients most commonly undergo conscious sedation, monitor patient vitals. The physician may also elect to apply a gel anesthetic in back of throat.

SET-UP

1. Install, power up and set up Mederi RF Generator as outlined in Operator's Manual provided with Generator.
2. Unpackage Stretta Catheter and Accessories.
3. Follow instructions in Operator's Manual for Catheter set up.

STRETITA PROCEDURE

The recommended treatments are 4 antegrade rings in the distal LES (8 treatment sites per ring) - levels 1-4, and 2 pull-back rings in the gastric cardia (12 treatment sites per ring) - levels 5-6.

1. Endoscopic inspection, measurement and Guidewire placement confirms the location and depth of the patients Z-Line (squamo columnar tissue) versus the fixed oral bite-block. The Z-line will serve as the reference landmark for the first 4 (Antegrade) of 6 total treatment levels.
2. From Standby mode, push the Power On/Mode button to advance to Ready mode, the Catheter is inserted to the Z-Line measurement, retracted 3 cm, and advanced to a position 1 cm above the Z-line for treatment level 1. Placement is confirmed versus the bite-block.
3. Attach suction to Catheter.
4. The PRV is attached to the air syringe. 30 ml of air is drawn into the syringe and inserted into the insufflation port on the Stretta Catheter until the PRV releases excess air. Inflate the balloon, then extend the needles.

ANTEGRADE TREATMENTS

5. The needles are extended to the full extent, and then retracted to the treatment depth. Impedance readings are monitored on the Generator screen. Optimal impedance readings below 200 indicate proper placement of needles.
6. Once proper placement is confirmed the foot pedal is depressed once, beginning a one minute treatment cycle.

BETWEEN TREATMENTS: Retract the needles and deflate the balloon. Pull the catheter back up to 25 cm, rotate 45 degrees to the right and then advance Catheter down guidewire to desired depth as measured against the bite block. If resistance is experienced while pulling the catheter back, release suction and then re-establish suction prior to next treatment.

NOTE: Use the shaft of the Catheter to rotate, NOT the handle.

7. Reinflate balloon using the PRV to prevent over-inflation. The foot pedal is depressed and a second set of lesions is created, establishing the first antegrade ring of 8 lesions. Three more rings are created in this manner: 0.5 cm above the Z-line, at the Z-line, and 0.5 cm below the Z-line.

Be sure to use the BETWEEN TREATMENT instructions above after each set of lesions is created in the Antegrade levels.

CARDIA PULL-BACK/RETROGRADE TREATMENTS

NOTE: For these treatment levels do not pull back higher than 2 cm above the Z-line, and do not advance into the stomach more than 2 cm below the measured Z-line.

8. Remove PRV and Guidewire. Do not use for cardia pull-back treatments.
9. For treatment level 5, advance the catheter into the fundus of the stomach and inflate the balloon to 25 ml, retract against the hiatus to fit snugly. Extend the needles and deliver the first treatment cycle in this level. At the completion of the first treatment cycle, retract the needles, advance catheter into the stomach, rotate 30 degrees to the right, extend needles and repeat, then repeat cycle again 30 degrees to the left of first treatment.
10. For treatment level 6, retract needles and advance the catheter again into the fundus of the stomach and deflate the balloon to 22 ml, then retract against the hiatus to fit snugly, above treatment level 5. Extend the needles and deliver the first treatment cycle in this level. Following the first treatment cycle, retract the needles, advance catheter into the stomach, rotate 30 degrees to the right, extend needles and repeat, then repeat cycle again 30 degrees to the left of first treatment.
11. Following completion of 6 levels of treatment, the balloon is deflated and the catheter is removed. An endoscopic inspection of the treatment area should confirm safe delivery and completion of Stretta Therapy.
12. Disconnect Stretta Catheter and dispose of device.
13. Follow shut down procedure outlined in Generator Operator's Manual.

STRETITA POST-PROCEDURE PATIENT INSTRUCTIONS

Provide patient with the following Post-Procedure instructions and specific contact person in event of adverse symptoms.

- Avoid instrumentation of esophagus for at least one month after treatment, specifically the passage of a nasogastric tube.
- Patient may be cautioned to refrain from use of oral NSAIDs for a period of 14-days after Stretta procedure
- Use liquid acetaminophen with or without codeine or equivalent preparation for post-treatment analgesia.
- Instruct patient to crush all medications or use liquid preparations of all medications, where possible, for at least one month after treatment.
- Modify diet after Stretta procedure; full liquids for 24 hours, soft diet for 2-weeks.
- Patient should continue previous anti-secretory regimen for 2-months after Stretta procedure.

Notify Stretta physician if any of the following symptoms occur:

Fever	Vomiting	Increasing abdominal discomfort
Chest pain	Bleeding	Difficulty swallowing
Nausea	Shortness of Breath	Tachycardia (fast pulse)

SYMBOLS/ABBREVIATIONS

USE BY

DO NOT REUSE

SEE INSTRUCTIONS FOR USE

REF

CATALOG NUMBER

STERILE EO

STERILIZED USING ETHYLENE OXIDE

CAUTION

IRRIGATION (FLUID IN)

SUCTION (FLUID OUT)

BALLOON
INFLATION
PORT

THUMB
CONTROL

BALLOON


NEEDLES

SUCTION CONNECTION

IRRIGATION CONNECTION

CONNECTING CABLE CONNECTION

NEEDLE POSITIONS

**NEEDLES
RETRACTED**
RIDGE ON THUMB
CONTROL ALIGNS
WITH REAR DOT


**FULL
NEEDLE
EXTENSION**

**RF DELIVERY
POSITION**

FULL NEEDLE EXTENSION
THUMB CONTROL IS PUSHED
ALL THE WAY FORWARD

RF DELIVERY
THUMB CONTROL
ALIGNS WITH
FRONT DOT

© Mederi Therapeutics Inc. 2014. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form, or by any means: electronic, mechanical, magnetic, optical or otherwise, without the prior written permission of Mederi Therapeutics Inc., 800 Connecticut Ave., Norwalk, CT USA.

Mederi Therapeutics Inc. reserves the right to change its products and services at any time to incorporate the latest technological developments. This guide is subject to change without notice. Although this guide has been prepared with every precaution to ensure accuracy, Mederi assumes no liability for any damages resulting from the application or use of this information.

Mederi® and Stretta® are trademarks of Mederi Therapeutics Inc.

U.S. Patents 6,440,128; 6,056,744; 6,846,312; 6,971,395; 6,712,814; 6,009,877; 6,613,047; 6,044,846; 6,749,607; 6,092,528; 7,462,179; 6,254,598; 6,673,070; 7,125,407; 6,197,022; and other patents pending.

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